

## COLONY ASSAY

### 1. Drug treatment after plating cells

- Seed MEFs on 24-wells plates.
- Trypsinize 6-well (10 cm<sup>2</sup>) with 0.5 ml TVP (5 min 37 °C).  
Add 0.5 ml complete + β + LIF and resuspend cells.  
Count cells → 4 - 8x10<sup>6</sup> cells per ml.
- Seed 500 cells per 2 cm<sup>2</sup> well:  
1 μl cells in 4 ml complete + β + LIF for 8 wells of a 24-wells plate.  
Use Eppendorf multipet to pipet 0.5 ml medium with cells in each well.
- Next day:
  - . (incubate cells with 40 μM O<sup>6</sup>-benzylG for 1 h at 37 °C if methylating or chloroethylating drugs are tested).
  - . incubate with drug for 1 h at 37 °C.
  - . wash cells with PBS.
  - . add per well 0.5 ml complete + β + LIF.
- Count colonies after 3 days (200-250 colonies in well without drug).

### 2. Drug treatment before plating cells

- Seed MEFs on 24-wells plates.
- Trypsinize 6-well (10 cm<sup>2</sup>) with 0.5 ml TVP (5 min 37 °C).  
Add 0.5 ml complete + β + LIF and resuspend cells.  
Count cells → 4-8 x 10<sup>6</sup> cells per ml.
- Pipet 100-150 μl cells in 8 ml complete + β + LIF.
- Use Eppendorf multipet to pipet 0.5 ml medium with cells in Nunc vials.
- Add 10 μl drug of appropriate concentration to each vial.
- Incubate at 37 °C for 1 h.
- Transfer 10 μl of treated cells to 24 wells containing 0.5 ml complete + β + LIF.
- Count colonies after 3 days (200-250 colonies in well without drug).

---

O<sup>6</sup>-benzylG (20 mM)

O <sup>6</sup> -benzyl-Guanine	4.82 mg
DMSO	1 ml

---