

DIFFERENTIATION/CULTURE OF EMBRYOID BODIES

Stage 1: 4 days, culture of ES cells on non-coated TC surface in depleted medium

- Seed two times $2.5 \cdot 10^6$ ES cells in 2x 6 cm TC dishes (non-gelatin-coated) in 2x 5 ml depleted medium (DM).
- Refresh medium every day.
- After 3 days of culture loosened parts of ES colonies which piled up the days before will float in medium; to prevent destruction and loss of these cell clumps refresh cells as follows:
 - . transfer medium with wide bore 10 ml pipette to 15 ml conical tube.
 - . let sediment ES cell clumps for 1 min.
 - . aspirate supernatant carefully, leave some medium.
 - . add fresh 5 ml DM, resuspend pellet carefully again with wide bore pipette and seed cell clumps in same dish.
- On this same day prepare agarose layers for stage 2 (next day):
 - . autoclave 2% agarose in PBS.
 - . wet surface of 2x 3 10 cm TC dishes with \pm 5 ml warm agarose solution each.
 - . remove air bubbles; leave 10' RT to get the thin layers.
 - . equilibrate each agarose layer ON (at least 4 hours) in the incubator with \pm 10 ml DM.

Stage 2: 6 days, culture on agarose layer in depleted medium

- Collect cell clumps of each of both dishes:
 - . aspirate sup with wide bore pipette, transfer to 15 ml conical tube.
 - . direct a jet of 5 ml DM across the surface of the dish to dislodge small cell clumps and transfer medium with loosened cell clumps to another conical tube; repeat this until the surface looks empty.
 - . let cell clumps sediment for 5' at RT.
 - . aspirate sup (carefully, pellet is very loosely packed).
- Aspirate medium of agarose coated dishes.
- Resuspend cell clumps in 20-30 ml DM and spread them over 2-3 agarose coated dishes: incubate untouched for 72 hrs.
- Grown cell clumps form embryoid bodies (EB) which do not attach to the agarose layer.
- After 3 days EB float in medium, refresh (aspirate medium plus EB, let sediment for 5' at RT, aspirate most of the sup, resuspend EB in \pm 10 ml DM per 10 cm dish and transfer to same agarose layer) and leave 3 days in incubator.
- In these days differentiated tissues can be seen in EB, after these three days (= after 10 days of culture) beating of heart tissue can be seen.

Stage 3: 4-5 weeks, culture on non-coated TC surface in complete medium

- Collect EB like before, but:
 - . Resuspend EB in CM.

. Depending on the number of EB, reseed on 2 - 4 x 6 cm TC dishes (not gelatin coated) with EB of 1 x 10 cm agarose coated petri dish (5 ml/6 cm petri dish).

- Leave dishes 2 days in incubator untouched to allow attachment of aggregates to the dishes (aggregates unattached after 2 days are unlikely to do so).
- Refresh medium of attached EB (thereby removing unattached ones).
- Change medium when necessary (approx. every 2 days) and observe the progress of differentiation over a period of 4-5 weeks.