

LACZ STAINING

1. In tissue culture

- Wash cells twice with PBS.
- Fix the cells 15 min at RT in PEMS + 4% formaldehyde.
PEMS:
 - 0,1 M PIPES pH 6.9
 - 2 mM MgSO₄
 - 1 mM EGTA
- Wash twice with PBS.
- Incubate at 37 °C in humidified incubator with LACZ stain.

LACZ stain:

- 50 mM NaPi pH7.3
- 150 mM NaCl
- 1.3 mM MgCl₂
- 3 mM K₄Fe(CN)₆
- 3 mM K₃Fe(CN)₆
- 1 mg/ml X-gal (20 mg/ml stock in dimethylformamide, at -20 °C)

NOTE: Reaction is light sensitive, therefore cover the tissue culture plates with aluminium foil.

2. *In situ* (Elles H3)

- Prepare mouse tissue on dry ice.
- Store at -80 °C if necessary.
- Cut cryosections on coated slides.

- Thaw slides ± 5 min.
- Allow fixation for 5 min. on ice in
 - 0.2% Glutaraldehyde
 - 0.1 M Phosphate buffer, pH 7.3
 - 5 mM EGTA
 - 2 mM MgCl₂
- Rinse twice in PBS + 2 mM MgCl₂ at R.T. (one rinse = 5 dips).
- 5 min in detergent rinse at R.T.
- Circle section with Dakopen.
- Apply on section 100-200 µl staining solution.

Staining sol.: - 20 mM Tris, pH 7.3
- 1 mg/ml X-gal

- staining sol. should be at R.T. before adding X-gal
- filtrate through 0.22 µm

- Incubate O.N. at 37 °C in humidified incubator.
- Rinse twice in PBS + 2 mM MgCl₂.
- Rinse once in aqua dest.
- Stain for 2 min in neutral red.
- Rinse three times in aqua dest.
- Rinse once in 70% ethanol.
- Rinse once in 100% ethanol.
- Air dry for 15 min.
- Rinse twice in xylene.
- Leave slides in xylene and mount one by one in DEPEX.

3. Whole mount

- Rinse embryos six times in 0.5 ml PBS.
- Fix embryos for 20 min (longer depending on size) in cold PBS containing:
 - 2% formaldehyde
 - 0.2% glutaraldehyde
 - 0.4% PVP-40
- Rinse embryos again extensively in PBS.
- Transfer embryos to 0.5 ml lacZ stain:
 - PBS with:
 - 5 mM K₄Fe(CN)₆
 - 5 mM K₃Fe(CN)₆
 - 25 mM MgCl₂
 - 2 mg/ml X-gal (Bethesda Res. Labs)
- Incubate in the dark for 16 to 24 h at 37 °C in humidified chamber.