

## LIPOFECTION

- Dissolve Tfx™ -10, -20 or -50 in 400  $\mu$ l nuclease-free water. Vortex and incubate at 65 °C for 1 min.

Store at -20 °C.

- Culture cells to 80% confluency on 24-well plate.
- Prepare complete medium without serum.  
Add DNA to medium, mix well  
Add Tfx™ at ratio Tfx™ : DNA = 3 : 1 according to table:

	Amount of DNA per well			
	0.25 $\mu$ g	0.50 $\mu$ g	0.75 $\mu$ g	1.00 $\mu$ g
medium	1.4 ml	1.4 ml	1.4 ml	1.4 ml
DNA	1.8 $\mu$ g	3.5 $\mu$ g	5.3 $\mu$ g	7.0 $\mu$ g
Tfx™	7.9 $\mu$ l	15.8 $\mu$ l	23.6 $\mu$ l	31.5 $\mu$ l

(ratio Tfx™ : DNA = 2 : 1 or 4 : 1 can also be tried)

- Incubate medium/DNA/Tfx™ lipofection mix at RT for 15 min.
- Incubate cells in 200  $\mu$ l lipofection mix at 37 °C for 30 min.
- Add 1 ml complete medium (with serum).
- Incubate cells for 24-48 h.
- A LacZ reporter plasmid can be used to measure transfection efficiency.

## NOTE

Optimal conditions for ES cells: 0.75-1.00  $\mu$ g of DNA per well  
Tfx™ -50 : DNA = 3:1